



Registered NDIS Provider


HANDICARE 1000 STAIRLIFT

FOR INDOOR &
OUTDOOR USE



HANDICARE 1000 STAIRLIFT

The Handicare 1000 stairlift has one of the **slimmest straight stairlift tracks on the market**. With no exposed gearing, the track provides a sleek look and takes up little room along your staircase, ensuring room for other stair users to get by comfortably and safely.





The freedom to live with independence in your own home.

- A standard weight limit of 137kg (200kg option also available)
- A slim profile of only 385mm (when seat, arms and footrest are folded up) saves space for other stair users
- Safety sensors immediately stop the lift if anything is detected in its path
- An adjustable seat height ensures the perfect fit for each user
- A manual swivel seat provides for safer transfer onto and off the lift (a powered upgrade option is available on the indoor model)
- A vinyl seat provides for easy cleaning
- A large, carpeted footrest gives greater traction
- Battery operated means the lift will still run in the event of a power failure (for a limited number of travels).



INCLUSIONS

- Retractable seatbelt
- Fold up seat, armrests and footrest
- Key operated switch (prevents unauthorised use)
- Manual emergency stop
- Soft start and stop
- Footrest and power pack safety edges
- Two remote controls

OPTIONAL EXTRAS

- 200kg limit
- Powered swivel seat *
- Powered footrest *
- Body harness
- Foot restraints
- Adjustable distance between armrests
- Wander lead (nurse control)
- Manual or powered hinge track *
- Perch Seat *

* For safety reasons these options are not available on the outdoor model.

SEAT OPTIONS

Two seat options are available to suit staircase designs and user requirements.

Smart Seat

The Smart Seat is our most versatile seat that offers users an array of options and upgrades:

- Easy to clean, durable and water resistant
- Adjustable seat height and width
- Easily upgraded to a power swivel seat

Perch Seat

Available on indoor model only.

- Ideal for those with knee or hip flexion problems, the Perch Seat alleviates the need to sit down each time the user rides the stairlift. As the seat is not as deep, but higher than a traditional seat, it allows the user to use the stairlift without adopting a fully seated position. Instead, the user is almost in a standing position, using the half seat for support
- The Perch Seat is also a perfect solution for narrow staircases where standard stairlift seats do not have clearance, or if the user's knees may touch the opposite banister
- Easy to clean, durable and water resistant
- A weight limit of 137kg

DESIGN FEATURES

Our stairlifts afford a variety of options with different levels of automation and features to meet your physical needs, personal preferences and budget. Choose from a range of options to create a stairlift that is right for you and your home.



Powered Swivel Seat

For an easier transfer, swivel your seat towards the landing once you reach the top (manual swivel is standard or upgrade to our automatic powered option).



Safety Sensors

The Handicare 1000 will stop if it meets an obstruction.



Easy-To-Use Toggle Control

Designed in collaboration with occupational therapists, the toggle is easy to use.



Perch Seat

The Perch Seat is ideal for users with limited hip and knee flexibility.



Hinged Track

A powered or manual hinged track is a smart solution to prevent the track from obstructing the doorway at the base of the stairs.






Waterproof Cover

The outdoor stairlift is UV protected and includes a waterproof cover so it is always dry and ready for use.





SEAT SPECIFICATIONS



	Smart Seat	Perch Seat
A Minimum swivel radius	670mm	
B Wall to inner seat back *	75mm	75mm
C Seat depth	360mm	320mm
With extended seat	425mm	
D Minimum folded width to footplate *	385mm	385mm
E Minimum open width to edge footplate *	650mm	650mm
F Armrest width – external	550mm	550mm
With armrest spacers	590mm	
With extended armrest spacers	635mm	
G Armrest width – internal	465mm	465mm
With armrest spacers	520mm	
With extended armrest spacers	565mm	
H Seat back height	455mm	455mm
I Footplate to seat height	465 - 540mm	675 - 875mm
J Minimum footplate height	70mm	70mm
Minimum track intrusion into staircase	115mm	115mm
Maximum weight limit - Standard version	137kg	137kg
Heavy duty version (XL) **	160kg	N/A
Extra heavy duty version (XXL) **	200kg	N/A
Speed	0.07m/sec	0.07m/sec
Staircase angle	27 - 56°	27 - 56°

* +45 mm for extra heavy duty (XXL) 2-way swivel with extended armrest spacers

** Not available with hinged rail

Colour Choices

Indoor - available in 6 leather grain, anti-slip upholstery options including;



Cocoa



Ruby



Sand



Sapphire



Slate



Jade

Outdoor - available in leather grain, anti-slip upholstery in 'Slate'.



Slate



How do I choose the right stairlift?

To ensure that your stairlift suits your needs, one of our team can visit your home to assess your staircase and discuss your requirements. This short consultation is without any obligation, and will be conducted by one of our qualified installers.



Each stairlift is designed to suit your individual staircase.

Using our specialist PhotoSurvey software we take measurements of your staircase, and then create an accurate 3D model – this software ensures the stairlift track is designed to perfectly fit your space. We also advise of any adjustments required to fit your stairlift, for example moving a handrail, and which side of the staircase would best suit your stairlift.



What are my stairlift options?

Stairlifts are either manufactured for straight or curved staircases, and for internal and external applications. Each of these stairlifts provide a range of features and functionality, and depending on your needs and budget there are options to increase the automation of the lift to make it even easier to use.

Visit our showroom today, call us on **1800 940 575** or visit **nationalstairlifts.com.au** for more information.

This publication provides general information only. No representation or warranties are made regarding the accuracy, completeness or relevance of the information provided. All measurements quoted are approximate values only and can vary. National Stairlifts accepts no liability for any loss or damage resulting from any reliance on the information provided in this publication. National Stairlifts is a registered business name.

© May 2019 National Stairlifts Pty Ltd ABN 88627357386. All rights reserved. This guide or any part of it may not be reproduced without prior written consent from National Stairlifts Pty Ltd. 'I support NDIS' and 'Registered NDIS Provider' trademarks owned by the Commonwealth Government.



GET IN TOUCH

T 1800 940 575
E info@nationalstairlifts.com.au
W nationalstairlifts.com.au

ABN 88627357386

National Stairlifts Pty Ltd

Unit 15
40 Counihan Rd
Rocks Riverside Industry Park
Seventeen Mile Rocks
QLD 4073