

Handicare Freecurve



1800 940 575
Unit 15/40 Counihan Road
Seventeen Mile Rocks QLD 4073

The Handicare Freecurve

Handicare has a solution for all situations and staircases.

All of our stairlifts are quiet, safe and easy-to-use. No matter what your staircase looks like, our compact and customised stairlifts can be installed anywhere. Even on the narrowest of staircases, there is enough walking space available after the stairlift is installed.

High quality in all aspects makes a Handicare stairlift a true asset to your independence. It

enables you to keep living in your own home, free from the risks associated with the stairs. A stairlift provides freedom of movement and peace of mind. Many of our satisfied customers experience that every day.




COMFORT

All of our stairlifts are fitted with a comfortable seat, armrests and a safety belt.



EASY-TO-USE


Ergonomic and intuitive to use, the joystick smoothly operates the stairlift.



CUSTOMISED


Each stairlift rail is designed and manufactured specifically to fit your staircase, ensuring a tight fit to maximize space for other stair users.

The Handicare Freecurve




COMPACT AND FOLDABLE

Fold up the stairlift when not in use to save space in corridors or hallways.



REMOTE CONTROL

Two remote controls are included for both upstairs and downstairs so you can send for your stairlift at the touch of a button.



ACTIVE SEAT

This option gently helps you get in and out of your stairlift, reducing strain on your joints and muscles.



SAFETY FIRST

Powered options eliminate the need to bend, safely folding or unfolding the footrest for you. The footrest is illuminated to prevent a tripping hazard in low-lit areas.

The advantages of your stairlift

With its many features, the Handicare Freecurve stairlift will be a valuable addition to your home.

Space

Once installed, there will still be enough free space on the stairs to walk up and down safely.



Safety stop

The stairlift is fitted with safety sensors that will stop the lift whenever it runs into an obstacle



Hinged rail




Dropnose



Compact rail endings

To avoid an obstruction at the bottom of the stairs (e.g. where there is a doorway), an automatic hinged rail with safety system or a dropnose is available.
Hinged rail Dropnose



Seat swivel

Swivel your stairlift seat toward the landing for a much easier and safer exit. This can be done manually or automatically. The optional Turn & Go enables Handicare to service even the narrowest stairs.



Active seat

Normally only available on lift chairs, gently supports your movements while sitting down or standing up.

A stairlift to co-ordinate with every interior

Of course it is important that your stairlift fits within the rest of your home décor. That is why we offer you the chance to choose from three seat styles with different varieties of upholstery and three standard rail colours.



CREAM WHITE



BROWN



GREY

Alliance Seat



BORDEAUX RED



BEIGE



COCOA BROWN



Elegance Seat



BORDEAUX RED



BEIGE



COCOA BROWN



Classic Seat



BORDEAUX RED



BEIGE



SLATE

