




HANDICARE 2000 STAIRLIFT

FOR INDOOR AND OUTDOOR USE



Registered NDIS Provider

The Handicare 2000

Handicare has a solution for all situations and staircases.



All of our stairlifts are quiet, safe and easy-to-use. No matter what your staircase looks like, our compact and customised stairlifts can be installed anywhere. Even on the narrowest of staircases, there is enough walking space available after the stairlift is installed.

High quality in all aspects makes a Handicare stairlift a true asset for your independence. It enables you to keep living in your own home, free from the risks associated with the stairs. A stairlift provides freedom of movement and peace of mind. Many of our satisfied customers experience that every day.


COMFORT

All of our stairlifts are fitted with a comfortable seat, armrests and a safety belt.



EASY-TO-USE

Ergonomic and intuitive to use, the joystick smoothly operates the stairlift.



CUSTOMISED

Each stairlift rail is designed and manufactured specifically to fit your staircase, ensuring a tight fit to maximise space for other stair users.



The Handicare 2000



COMPACT AND FOLDABLE

Fold up the stairlift when not in use to save space in corridors or hallways.



REMOTE CONTROL

Two remote controls are included for both upstairs and downstairs so you can send for your stairlift at the touch of a button.



RAIL START AND END

The rail start and ending can be adjusted to your personal preference and situation. A flush finish, for example, is ideal if there isn't a lot of space at the top of the stairs.



COMFORT AND SAFETY

The powered footrest eliminates the need to bend, safely folding or unfolding the footrest for you.

The advantages of your stairlift

With its many features, the Handicare 2000 stairlift will be a valuable addition to your home.

Space

Once installed, there will still be enough free space on the stairs to walk up and down safely.



Safety stop

The stairlift is fitted with safety sensors that will stop the lift whenever it runs into an obstacle.



Automatic hinged rail

If there is limited space at the bottom of the stairs or if the rail blocks the hallway or a door, you can choose the automatic hinged rail.

Seat swivel

Upon arriving at the top of the stairs, you can swivel the seat towards the top landing for a much easier and safer exit. This can be done manually but automatic seat swivelling is available as an option as well. Simply hold the ergonomically designed joystick towards the top landing and the seat will swivel automatically.



Load capacity

The weight limit on the Handicare 2000 is 115 kg. The weight limit can be increased to 137 kg with the heavy duty kit.

A stairlift in a block of units?


The Handicare 2000 is the right choice.

Handicare always strives to offer the right stairlift for every possible situation. This is one of the reasons why we are a leader in the installation of stairlifts in units and have a large number of satisfied customers.

- Our developers have given thought to the specific requirements of stairlifts in units and have placed special emphasis on comfortable, safe and functional solutions.
- The stairlift is supplied with key switch and key as standard to prevent unauthorised use.
- It is possible to stop the stairlift on each floor with additional parking positions.
- Easily call the stairlift or send it to the next floor with the remote controls.
- In the event of a power failure the stairlift can still be used for a certain number of trips thanks to built-in batteries.



National Stairlifts will gladly assist you with the necessary steps for installing a stairlift in a unit.



A stairlift to co-ordinate with every interior


Of course it is important that your stairlift complements the rest of your home décor. That is why we offer you the chance to choose from four seat styles with different types of upholstery and three standard rail colours.



GUN METAL



BROWN



LIGHT GREY

Smart Seat



RUBY



CACAO



SAND



SAPPHIRE



SLATE



JADE



Simplicity Seat



SAND

Style Seat



IVORY
WHITE



CHERRY
RED



VANILLA
BEIGE



Smart Perch Seat



RUBY



CACAO



SAND



SAPPHIRE



SLATE



JADE



If you have a narrow staircase or problems with your knee or hip joints and cannot sit down, the Smart Perch seat is an ideal choice for you. The raised seat and the reduced seat depth allow you to use the lift in a standing position while being secured by a seatbelt.

How do I choose the right stairlift?

To ensure that your stairlift suits your needs, one of our team can visit your home to assess your staircase and discuss your requirements. This short consultation is without any obligation, and will be conducted by one of our qualified installers.



Each stairlift is designed to suit your individual staircase.

Using our specialist PhotoSurvey software we take measurements of your staircase, and then create an accurate 3D model – this software ensures the stairlift track is designed to perfectly fit your space. We also advise of any adjustments required to fit your stairlift, for example moving a heater, and which side of the staircase would best suit your stairlift.



What are my stairlift options?

Stairlifts are either manufactured for straight or curved staircases, and for internal and external applications. Each of these stairlifts provide a range of features and functionality, and depending on your needs and budget there are options to increase the automation of the lift to make it even easier to use.

Our new mobile showroom enables clients to see first-hand how easy our stairlifts are to use, and how life changing they can be.

We provide 4D imagery of the stairlift in their house, and an accurate quotation on supply and installation.



Please contact us today to request further information on a suitable stairlift for your home. We provide free, no-obligation quotations.

Visit our Brisbane showroom or request a visit from our mobile showroom
Call 1800 940 575
info@nationalstairlifts.com.au
nationalstairlifts.com.au

This publication provides general information only. No representation or warranties are made regarding the accuracy, completeness or relevance of the information provided. All measurements quoted are approximate values only and can vary. National Stairlifts accepts no liability for any loss or damage resulting from any reliance on the information provided in this publication. National Stairlifts is a registered business name.

© September 2020 National Stairlifts Pty Ltd ABN 88627357386. All rights reserved. This guide or any part of it may not be reproduced without prior written consent from National Stairlifts Pty Ltd. 'I support NDIS' and 'Registered NDIS Provider' trademarks owned by the Commonwealth Government.

QBCC 15102469

QLD Electrical Licence 85362
NSW Electrical Licence 331996C
VIC Electrical Licence 28727
ACT Electrical Licence 20181244
NT Electrical Licence C3911



GET IN TOUCH

T 1800 940 575
E info@nationalstairlifts.com.au
W nationalstairlifts.com.au

ABN 88627357386

National Stairlifts Pty Ltd

15/40 Counihan Rd
Rocks Riverside Industry Park
Seventeen Mile Rocks
QLD 4073